DEVELOPMENT AND APPLICATION OF REAL-TIME PCR METHODS FOR THE DETECTION AND TYPING OF EXTENDED-SPECTRUM **b**-LACTAMASES IN CLINICAL ISOLATES OF ENTEROBACTERIACEAE

*M.V. Edelstein, A.N. Ekimov, O.V. Shevchenko Inst. of Antimicrobial Chemotherapy, Smolensk, Russia* 

#### РАЗРАБОТКА И ПРИМЕНЕНИЕ МЕТОДОВ ПЦР В РЕЖИМЕ РЕАЛЬНОГО ВРЕМЕНИ ДЛЯ ИДЕНТИФИКАЦИИ Ъ-ЛАКТАМАЗ РАСШИРЕННОГО СПЕКТРА У КЛИНИЧЕСКИХ ШТАММОВ ЭНТЕРОБАКТЕРИЙ

М.В. Эйдельштейн, А.Н. Екимов, О.В. Шевченко НИИ антимикробной химиотерапии, Смоленск, Россия



## EXTENDED-SPECTRUM b-LACTAMASES (ESBLs)

- Active against all b-lactams (penicillins, I-IV gen. cephalosporins, aztreonam) except cephamycins and carbapenems
- (-<sup>-</sup>) Inhibited by active-site directed inhibitors (clavulanic acid, sulbactam, tazobactam)
- Ambler molecular class A or class D (active site serine b-lactamases)
- ➢ Bush, Jacoby and Medeiros functional groups: 2be, 2d, 2f
- > Multiple genetic types: CTX-M, SHV, TEM, etc.
- > Typically plasmid-encoded, often associated with mobile elements
- Produced mainly by Enterobacteriaceae

## WHY ESBLS ARE IMPORTANT?


- Confer resistance to the most commonly used antibiotics which constitute the first-line therapy for nosocomial infections
- Difficult to detect by routine susceptibility tests:
   FALSE susceptible ESBL producers are a common cause of treatment failures with modern cephalosporins
- Complex epidemiology: may rapidly spread by clonal transmission or plasmid transfer
- Frequently associated with resistance determinants to non-blactam agents (fluoroquinolones, aminoglycosides, tetracyclines, sulfonamides, etc).
  ESBL producers are often multiply resistant

### **MOLECULAR TYPING OF ESBLs**


- Provides a basis for understanding the epidemiology and evolution of ESBLs
- There is no one single molecular test to detect all ESBL types
- Numerous methods have been proposed for rapid identification of ESBLs that belong to certain genetic groups (e.g. TEM, SHV, CTX-M):
  - Spoligotyping; PCR-RFLP;
  - PCR-SSCP; LCR;
  - minisequencing; real-time PCR
- But none of these techniques can identify all members of any group in a single step

#### **THE SHV-TYPE ESBLs**

- SHV-1 penicillinase (»32 kDa) is a direct progenitor
  - species-specific (chromosomally encoded) in Klebsiella pneumoniae
  - often plasmid mediated in different species of *Enterobacteriaceae*
- > SHV ESBLs evolve by acquisition of point mutations
- Over 60 amino acid variants currently exist\*




## **THE SHV-TYPE ESBLs**




Substitutions at only a few amino acid residues are required for extended-spectrum activity

#### MGB ECLIPSE<sup>™</sup> PROBES\*




- Contain minor grove binder (MGB) and Eclipse dark quencher at 5' end
- Contain fluorescent dye at 3' end
- Short probes(9-15 nt) with high Tm and precise binding to target


#### MGB ECLIPSE<sup>™</sup> PROBES\*



- Quenched in random coil state
- Fluorescing when bond to target
- 5'-nuclease resistant **P** suitable for postamplification meting curve analysis



## DESIGN OF MULTIPLEX SINGLE-TUBE PCR AND MELTING CURVE ANALYSIS FOR DETECTION OF SHV ESBLs




- > 4 differentially labelled probes complementary to WT sequences at key mutation sites
- Each mutation specifically shifts the Tm of the corresponding probe
- Asymmetric primer ratio to favour formation of the strands complementary to the probes

A.Ekimov et al., 14th ECCMID, 2004, P944


## **CONTROL STRAINS WITH KNOWN MUTATIONS IN SHVs**

Mutation	Detected with probe
WT, non-ESBL control (SHV-1)	All probes
Ser-238 (SHV-2, -3, -4)	238-FAM
Ser-238 + Lys-240 (SHV-5)	238-FAM
Ala-238 (SHV-18)	238-FAM
Ala-179 (SHV-6)	179-TET
Asn-179 (SHV-8)	179-TET
Gly-179 (site-directed mut.)	179-TET
Asp-156 (site-directed mut.)	156-Cy3.5
Ser-149 (site-directed mut.)	146-149- <b>Cy</b> 5
Val-146 (site-directed mut.)	146-149- <b>Cy</b> 5

## DETECTION AND DIFFERENTIATION OF ALL THE KNOWN ESBL MUTATIONS BY MELTING CURVE ANALYSIS




## **ANALYSIS OF CLINICAL ISOLATES**




Rotor-Gene 2000 Real-Time PCR System

## **ANALYSIS OF CLINICAL ISOLATES**



Rotor-Gene 2000 Real-Time PCR System

## **ANALYSIS OF CLINICAL ISOLATES**




Rotor-Gene 2000 Real-Time PCR System

## DETECTION OF SHV ESBLs BY MULTIPLEX REAL-TIME PCR WITH MGB ECLIPSE™ PROBES


- All the known mutations conferring extended-spectrum activity on SHV b-lactamases can be detected and discriminated in a single reaction
- To our knowledge, this is the first example of detection of 8 mutations in 5 codons using a single-tube real-time PCR
- New mutations at the key codons can also be identified
- Plasmid-mediated SHV ESBLs can be detected on the background of chromosomally-encoded SHV-1 in K.pneumoniae
- Extremely fast and processive
   (< 4 h for DNA isolation, PCR and melting curve analysis with</li>
   33 strains without opening the tube)
- > No risk of contamination by PCR products

## THE CTX-M-TYPE ESBLs



http://www.lahey.org/studies/other.asp


## REAL-TIME PCR WITH UNIVERSAL CTX-M PRIMERS AND SYBR GREEN I




- Primers were designed to amplify 519bp PCR-product from all the subtypes
- Members of different clusters can be distinguished by melting curve analysis

M.Edelstein et al., 15<sup>th</sup> ECCMID, 2005, P. 649

#### INCIDENCE OF ESBL-PRODUCING ENTEROBACTERIACEAE IN RUSSIAN ICUs



M.Edelstein et al., ICAAC, 2004, P.C2-1331




\* S.K. Bouchillon et al., IJAA 2004 (24): 119–24

\*\*M.Edelstein et al., ICAAC, 2004, P.C2-1331

## DRAMATIC INCREASE IN THE PROPORTION OF CTX-M ESBLs IN RUSSIA

**1998-99 – NPRS study\*** 


2003 – RESORT study\*\*



\* *M.Edelstein et al., ICAAC, 2004, P.C2-1331* 


# SPECIES DISTRIBUTION AMONG ESBL-PRODUCING ENTEROBACTERIACEAE

**RESORT** study – 2003 **ESBL-producing** strains (n=384)



M.Edelstein et al., ICAAC, 2004, P.C2-1331

## PREVALENCES OF DIFFERENT ESBL TYPES AND THEIR COMBINATIONS IN RUSSIAN NOSOCOMIAL ISOLATES



## **OUR THANKS TO:**

Dr. Antone Medeiros (Miriam Hospital, Providence, RI, USA)

Dr. Marek Gniadkowski (National Institute of Public Health, Warsaw, Poland)

Dr. Rafael Canton (Hospital Universitario Ramon y Cajal, Madrid, Spain)

Prof. Herbert Hächler (Institute of Veterinary Bacteriology, University of Bern, Switzerland)

> Prof. Adolf Bauernfein (MICOER Institute, Munich, Germany)

for providing reference **b**-lactamase producing strains

Dr. Evgeniy Belousov (Epoch Biosciences, Bothell, WA, USA) for assistance in designing the MGB Eclipse probes

This work was supported by the Russian Foundation for Basic Research and Administration of Smolensk Region